

SURGICAL DRAPES  
TEST OF TEAR STRENGTH OF SURGICAL DRAPES

TEST STANDS (MECHANIC) | KERN TVL

Test of tear strength of surgical drapes. According to the regulations of EN 13795, surgical drapes must have a high strength (wet and dry). The tear resistance test can be used to check whether the fabric meets these requirements



Model	Max. Force	Max. Travel distance	Price
TVL		230 mm	€ 450,-



<https://partner.kern-sohn.com/Segment/B0253>

\* DAkkS calibrations are also possible on site  
All prices are net, valid until 20.01.2026 plus shipping cost

Personal consultation from our product specialist:



**Irmi Russo** Productspecialist  
Tel : +49 7433 9933-208  
[russo@kern-sohn.com](mailto:russo@kern-sohn.com)

Application number: EN - 6 13

**KERN & SOHN GmbH**  
Ziegelei 1  
72336 Balingen, Germany  
[info@kern-sohn.com](mailto:info@kern-sohn.com)  
[www.kern-sohn.com](http://www.kern-sohn.com)

