

UTILITY CERAMICS MANUFACTURER

PRECISE WEIGHING OF INDIVIDUAL DOSES OF THE GLASS GRANULATE, CALLED "GLASS FRITS"

PALLET SCALES | KERN VHB

Precise weighing of individual doses of the glass granulate (glass frits) from different silos in portions of 50 kg to 150 kg. These precisely weighed amounts are used to mix a suspension for the enamelling process



Model	Weighing capacity	Readability	Net weight	Price
VHB 2T1	2000 kg	1 kg	106 kg	€ 1,930.-
<b>Suitable Accessories:</b>				
● DAkkS calibration certificate 963-131 - 963-131 - €290.-				



<https://partner.kern-sohn.com/Segment/B092E>

\* DAkkS calibrations are also possible on site  
All prices are net, valid until 03.11.2025 plus shipping cost

Personal consultation from our product specialist:



**Jürgen Fünfer**  
Productspecialist  
Tel : +49 7433 9933-245  
juergen.fuenfer@kern-sohn.com

Application number: EN - 614

**KERN & SOHN GmbH**  
Ziegelei 1  
72336 Balingen, Germany  
info@kern-sohn.com  
www.kern-sohn.com