



## TENNIS SPORTING GOODS RETAIL

### TESTING AND RECALIBRATION OF A STRINGING MACHINE FOR TENNIS RACKETS

#### DIGITAL FORCE-MEASURING DEVICES (INTERNAL SENSOR) | KERN FK

The stringing of the racket is an important criterion for tennis players and has a significant influence on the playing characteristics. Every tennis player not only improves his game with the right covering, but also protects his health. That is why the right clothing is extremely important and the correct pulling power of the clothing machine is to be checked regularly



Model	Max. Force	Readability	Price
FK 50	50 N	0.02 N	€ 275.-
<b>Suitable Accessories:</b> <ul style="list-style-type: none"><li>• DAkkS calibration certificate 963-2610 (force direction compression   measurement result unit N   calibration standard DKD-R 3-3 - sequence A) - 963-2610 - €345.-</li></ul>			



<https://partner.kern-sohn.com/Segment/B0B15>

\* DAkkS calibrations are also possible on site  
All prices are net, valid until 19.11.2025 plus shipping cost

Personal consultation from our product specialist:



**Irmi Russo** Productspecialist  
Tel : +49 7433 9933-208  
[russo@kern-sohn.com](mailto:russo@kern-sohn.com)

Application number: EN - 658

**KERN & SOHN GmbH**  
Ziegelei 1  
72336 Balingen, Germany  
[info@kern-sohn.com](mailto:info@kern-sohn.com)  
[www.kern-sohn.com](http://www.kern-sohn.com)