

SAVINGS BANKS

AUTHENTICATION OF GOLD BARS AND GOLD PIECES

MATERIAL THICKNESS GAUGES | KERN TN

Whether for pawn shops or for buying or selling, the authenticity of gold bars and pieces of gold is always of central importance. The authenticity can easily be determined with an ultrasonic measuring device and the appropriate software and protects against high financial risks.



Model	Max. Pulse-Echo	Price
TN GOLD 80	80 mm	€ 770.-
Suitable Accessories: <ul style="list-style-type: none">• Factory calibration certificate 961-113 - 961-113 - €167.-		



<https://partner.kern-sohn.com/Segment/B10CC>

* DAkkS calibrations are also possible on site
All prices are net, valid until 24.12.2025 plus shipping cost

Personal consultation from our product specialist:



Irmi Russo Productspecialist
Tel : +49 7433 9933-208
russo@kern-sohn.com

Application number: EN - 7 12

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com