



FARMER FOR MALTING BARLEY (TEST SEGMENT)
DETERMINATION OF THE RESIDUAL MOISTURE CONTENT OF PRE-DRIED BARLEY

MOISTURE ANALYSER | KERN DAB

After the barley has dried, samples are taken to determine the residual moisture content to protect against rapid germination or mold formation)



Model	Weighing capacity	Readability	Price
DAB 100-3	110 g	0.001 g	€ 1,250.-
Suitable Accessories: <ul style="list-style-type: none">• DAkkS calibration certificate 963-127 - 963-127 - €103.-			



<https://partner.kern-sohn.com/Segment/B000E>

* DAkkS calibrations are also possible on site
All prices are net, valid until 06.11.2025 plus shipping cost

Personal consultation from our product specialist:



Aleksandar Delic
Productspecialist
Tel : +49 7433 9933-120
aleksandar.delic@kern-sohn.com

Application number: EN - 89

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com