

OCCUPATIONAL THERAPISTS (DOCTORS FOR THERAPY BY SPECIFIC MOVEMENTS) MUSCLE STRENGTH TEST OF THE HANDS IN OCCUPATIONAL THERAPY

HAND GRIP DYNAMOMETERS | KERN MAP

The individual hand strength can be an indication of the health to be expected in the future. The hand strength can be checked efficiently with a force gauge - a quick first check of the patient.



Model	Weighing capacity	Readability	Price
MAP 80K1	80 kg	100 g	€ 285,-



<https://partner.kern-sohn.com/Segment/B100C>

* DAkkS calibrations are also possible on site
All prices plus VAT and shipping costs valid until 02.02.2026

For expert advice, please contact:



Irmi Russo Productspecialist
Tel : +49 7433 9933-208
russo@kern-sohn.com

Application number: EN - 956

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com