


ANIMAL SUPPLEMENTS MANUFACTURER
CONTROLLING WEIGHING OF CONTAINERS

PLATFORM SCALES | KERN IOC

Dietary supplements for animals in the form of pellets are packed by hand in 500g cans or 5000g plastic bags, weighed and prepared for dispatch. Exact control of the filling weight is important, since filling the container too low would result in complaints, and filling too high would reduce the return



Model	Weighing capacity	Readability	Price
IOC 6K-4	3 kg; 6 kg	0.0001 kg; 0.0002 kg	€ 390.-



<https://partner.kern-sohn.com/Segment/S10475>

* DAKS calibrations are also possible on site
All prices are net, valid until 05.11.2025 plus shipping cost

Personal consultation from our product specialist:



Jürgen Fünfer
Productspecialist
Tel : +49 7433 9933-245
juergen.fuenfer@kern-sohn.com

Application number: EN - 965

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen, Germany
info@kern-sohn.com
www.kern-sohn.com